





ADVANCE •

ADVANCE SERIES

Medical Simulators
Tailored to your needs



• ADVANCE TRAINING STATION



PORTABLE •

ADVANCE

Laparoscopic Simulator - LAPARO Advance

LAPARO Advance is the highest class of personal laparoscopy trainers. The device is portable, yields an incredible build durability and design quality. The full HD camera which can be rotated and adjusted allows to practice at all configurations.

ADVANTAGES

- Full HD 30° camera with adjustable ball joint
- 8 x Ø5 and 2 x Ø10 tool slots
- LED lighting for optimal illumination of the operating field
- Compatible with all LAPARO training modules
- Allows to track and record training
- A special solid casing designed for many hours of training
- Convenient in transport
- No subscription - buy once and practice as much as you want

ALL ADVANCE SETS INCLUDE

- LAPARO Advance simulator with build-in HD 30° and LED lighting
- USB with LAPARO software, compatible with Windows, MacOS and Linux
- Suturing Kit - Needle and a Thread x 2
- User Manual
- 2 year Worldwide Guarantee



ADVANCE MODULES

Laparoscopic Simulator - Advance Modules

- 6 Training Modules
 - Universal Holder
 - Marbles
 - Threading
 - Rubber Bands
 - Shapes/Suturing
 - Cutting/Suturing

EAN	Dimensions	Weight	Guarantee
5904024551010	22 x 33 x 45 cm	3,6 kg	2 years



ADVANCE COMPLETE

Laparoscopic Simulator - Advance Complete

- 4 training-purpose laparoscopic instruments
 - Dissector
 - Grasper
 - Scissors
 - Needle Holder
- 6 Training Modules
 - Universal Holder
 - Marbles
 - Threading
 - Rubber Bands
 - Shapes/Suturing
 - Cutting/Suturing

EAN	Dimensions	Weight	Guarantee
5904024551027	22 x 33 x 45 cm	4,4 kg	2 years



ADVANCE PEDIATRIC

Laparoscopic Simulator - Advance Pediatric

- 4 Pediatric Training laparoscopic Instruments
 - Dissector Ø3 mm
 - Grasper Ø3 mm
 - Scissors Ø3 mm
 - Needle Holder Ø3 mm
- 8 Training Modules
 - Universal Holder
 - Marbles
 - Threading
 - Rubber Bands
 - Shapes/Suturing
 - Cutting/Suturing
 - Holder with Adjustable Angle
 - Platform
- Trocars Ø3 mm x 2

EAN	Dimensions	Weight	Guarantee
5904024551034	22 x 33 x 45 cm	4,8 kg	2 years

